

Broadley
James®

BIO
NET

Bioreactors

You
Now
Have A
Choice!

toll-free: 800 . 288 . 2833

phone: 949 . 829 . 5555

fax: 949 . 829 . 5560

email: sales@broadleyjames.com


19 Thomas

Irvine, California 92618

www.broadleyjames.com



- Suitable for Cell Culture & Microbial Applications
- Like-for-Like Replacement of the Other Leading Manufacturer
- Single-Walled & Jacketed
- pH, DO, Temperature & Optical Density Sensors Available
- Direct Drive & Magnetically Coupled Agitation Assemblies
- Individual Bioreactor Components, Complete Assemblies, Rebuild Kits & O-Ring Packs Available
- Large Inventory Available for Immediate Delivery



FEATURES & BENEFITS

**Complete Bioreactor Systems in a variety of volumes:
3L, 5L, 7L, 15L & 20L.**

Single-Walled & Jacketed Styles

The BioNet® Bioreactor is the latest in autoclavable benchtop glass vessels. Available in a range of sizes: 3, 5, 7, 15 and 20 liters total volume. All metal components are manufactured from 316L stainless steel with electropolish.

Integral scaling, improved component designs and extensive product inventory makes this the most user-friendly bioreactor available.

Special configurations are available upon request.

Agitation motors, sensors, cables and control systems are available.

Dimensional and Data Sheets available at www.broadleyjames.com

For information about our Bioreactor Control System:

Please contact us at: **800.288.2833**